


Table of Contents

Entering in CCF/Specimen ID along with Cup Lot/Expiration	4
Donor Signature	6
Collectors Signature with Pin	7
Enter Results Page	8
Negative Results Page.....	9
Non-Negative Results Page.....	10
Submit a Collection Event Window (Non-Negative).....	11
Printing CCF1 Testing Facility Copy Window	12
CCF1 Laboratory Copy Example	13



Navigation

- Account Administration
- Occupational Health Screening
 - Schedule New Order
 - Pending Orders
 - New i3ccf Order**
 - i3ccf Orders
 - Completed Collections
 - Service Referrals
- Participants
- Help and Tutorials
- i3ccf compliance

Details

No details to display

New i3ccf Order

Use this wizard to order, collect, and result an i3ccf order. **All required fields are in bold**

Collection Information

Company:
ABC Company Demo Account

Company Location:
ABC Company Demo Account: ABC Company Demo Account

Package: 10 PANEL I3CCF **Reason for Test:** PRE-EMPLOYMENT

Order Expiration: 2018-01-19 23:59:59 Location Code: Federal Agency: Observed Collection Required:

Collection Site

Collection Site:
ABC Company Demo Account Collection Site -i3CCF Collection - 9501 Northfield Boulevard, , Denver, CO, 80238

Donor Information

First Name: John **Middle Name:** **Last Name:** Doe

SSN/EID: 55555555 **DOB:** **Email:**

Daytime Phone: 555-555-5555 **Evening Phone:**

Next >

New i3ccf Order

Use this wizard to order, collect, and result an i3ccf order. **All required fields are in bold**

Order Selection

Service: POCT Collection
Status: Ordered, Not Scheduled
Order Expiration: 2018-01-19 23:59:59
Reason For Test: PRE-EMPLOYMENT
Federal Agency:
Description: NexScreen NEXSCREEN CLIA 10 PANEL
Panel Codes: C-61012N
Laboratory: NexScreen
Collection Site: ABC Company Demo Account Collection Site -i3CCF Collection

Test Performing

NEXSCREEN CLIA 10 PANEL: Amphetamines, Barbiturates, Benzodiazepines, Cocaine, Methamphetamines, Methadone, Opiates, Oxycodones, Marijuana, MDMA

Employer

Location: ABC Company Demo Account
Address: 9501 Northfield Boulevard, Denver, CO 80238
Phone: 555-555-5555
Fax:

Donor

First Name: John	Middle Name:	Last Name: Doe
SSN/EID: 55555555	DOB:	Email:
Daytime Phone: 555-555-5555	Evening Phone:	

< Previous Save Save and Continue >

CONFIRM ORDER

ⓘ ATTENTION YOU ARE ABOUT TO CREATE A NEW ORDER.

Are you sure you wish to proceed?

Yes

No

Specimen Collection - John Doe - ***5555

Specimen Collection

Employer Information

Location: ABC Company Demo Account
Address: 9501 Northfield Boulevard, Denver, CO 80238
Phone: 555-555-5555
Fax:

MRO Information

MRO Name: David Nahin M.D.
Address: 9501 Northfield Blvd.Denver, CO 80238
Phone: 877-585-7366
Fax: 855-253-5666

Donor Information

First Name: John	Middle Name:	Last Name: Doe
SSN/EID: 55555555	DOB:	
Daytime Phone: 555-555-5555	Evening Phone:	

Verify Donor

Photo ID Employee Representative

Summary

Service: POCT Collection
Status: Ordered, Not Scheduled
Order Expiration: 2018-01-19 23:59:59
Reason For Test: PRE-EMPLOYMENT
Federal Agency:
Description: NEXSCREEN CLIA 10 PANEL
Panel Codes: C-61012N
Laboratory: NexScreen
Collection Site: ABC Company Demo Account Collection Site -i3CCF Collection

Next >

Entering in CCF/Specimen ID along with Cup Lot/Expiration

Specimen Collection - John Doe - *5555**

Collector and Site

Collector Name: KAREN GARCIA **Collector Phone:** 555-555-5555

Collection Site: ABC Company Demo Account Collection Site -I3CCF Collection - 9501 Northfield Boulevard, , Denver, CO, 80238

Instant Device

Enter the Specimen ID from the i3ccf Donor Consent & Label Form and enter the Lot Number and Expiration Date located on the device packaging.





CCF/Specimen ID: _____ **Lot:** _____

Product Expiration Date: _____

[Product Insert](#)
[Detailed Collection Steps](#)

Instant Cup Urine Collection

1. Give the Instant Test Cup device to the Donor
2. Have Donor void specimen in the Device in a secure bathroom facility

Suspend Collection < Previous Next >

Specimen Collection - John Doe - ***5555 - CCF/Specimen ID Number IS1234567890

Specimen Collection

Specimen Integrity and Validity

Did the donor **refuse** to comply with any part of the collection process? Yes No

Did the donor provide **sufficient quantity** for this specimen type? Yes No

Is the **temperature within acceptable range** (For urine specimens only; 90-100 degrees Fahrenheit)? Yes No Not Applicable

Was this collection an **observed** collection? Yes No




Is there any sign of specimen **tampering**? Yes No

Enter Remarks/Comments (**Notice: These remarks will be displayed on the CCF**):

Specimen Collection:

Single

Notes

 Add
 Search  Clear

Author	Date/Time	Subject	Type

Donor Signature

Specimen Collection - John Doe - ***5555 - CCF/Specimen ID Number IS1234567890

Specimen Collection

Employer Information Location: ABC Company Demo Account Address: 9501 Northfield Boulevard, Denver, CO 80238 Phone: 555-555-5555 Fax:	MRO Information MRO Name: David Nahin M.D. Address: 9501 Northfield Blvd.Denver, CO 80238 Phone: 877-585-7366 Fax: 855-253-5666
--	--

Collection Summary

Reason For Test: PRE-EMPLOYMENT
Federal Agency:
Drug Test to be Performed: C-61012N NEXSCREEN CLIA 10 PANEL
Laboratory:
Temperature In Range (Between 90-100): Yes
Collection: Single
Observed: No

Donor Signature

Donor Refused to Sign

I certify that I have provided the specimen(s); I have not adulterated them in any manner; each was sealed in my presence; and the information provided on this form and label(s) are correct.

Suspend Collection < Previous Next >

Collectors Signature with Pin

Specimen Collection - John Doe - ***5555 - CCF/Specimen ID Number IS1234567890

Collection Summary

Reason For Test: PRE-EMPLOYMENT
Federal Agency:
Drug Test to be Performed: C-61012N NEXSCREEN CLIA 10 PANEL
Laboratory:
Temperature In Range (Between 90-100): Yes
Collection: Single
Observed: No

Collector Signature

I certify that the specimen(s) given to me by the donor/participant identified in this form were collected, handled, and witnessed by me.

Click and drag using the left mouse button to sign.



Clear Use Topaz Pad

Collector Pin:


Suspend Collection < Previous Next >

Enter Results Page

Specimen Collection - John Doe - ***5555 - CCF/Specimen ID Number IS1234567890

Specimen Collection

Results Instructions



[Detailed Collection Steps](#)

[Product Insert](#)

Results Entry

Final Result: **Collector Pin:**

Enter Remarks/Comments (**Notice: These remarks will be displayed on the CCF**):

Delivery Service

Specimen Released To:


Notes

Negative Results Page

Specimen Collection - John Doe - ***5555 - CCF/Specimen ID Number IS1234567890

Specimen Collection

Results Instructions



[Detailed Collection Steps](#)
[Product Insert](#)

Results Entry

Final Result: Negative **Collector Pin:** ●●●●

Enter Remarks/Comments (Notice: These remarks will be displayed on the CCF):

Test Results

Drug Name	Result
Amphetamines	NEGATIVE
Barbiturates	NEGATIVE
Benzodiazepines	NEGATIVE
Cocaine	NEGATIVE

Suspend Collection
< Previous
Next >
View i3CCF2 - Donor/Customer Copy

SUBMIT A COLLECTION EVENT X

! **ATTENTION YOU ARE ABOUT TO SUBMIT A COLLECTION EVENT**

Donor Name: Doe, John

Date of Collection: Thu Jan 11 2018 09:47:52 GMT-0700 (Mountain Standard Time)

Collector Name: KAREN GARCIA

Product Name: POCT Collection

Description: NEXSCREEN CLIA 10 PANEL

RESULT = NEGATIVE

Confirm
Go Back

Result

Non-Negative Results Page

Specimen Collection - John Doe - ****5555 - CCF/Specimen ID Number IS123456789

Results Entry

Final Result: Non-Negative **Collector Pin:** ••••

Enter Remarks/Comments (**Notice: These remarks will be displayed on the CCF**):

Test Results

Drug Name	Result
Amphetamines	NEGATIVE
Barbiturates	NEGATIVE
Benzodiazepines	NEGATIVE
Cocaine	NEGATIVE
Marijuana	NEGATIVE
MDMA	NEGATIVE
Methadone	NON-NEGATIVE
Methamphetamines	NEGATIVE
Opiates	NEGATIVE
Oxycodones	NEGATIVE

Validity

Validity Name	Result

POCT Non-Negative Results Handling Procedures

All non-negative results are presumptive and should be confirmed by an alternate method. The preferred confirmation method is

Suspend Collection
< Previous
Next >
View i3CCF2 - Donor/Customer Copy

Submit a Collection Event Window (Non-Negative)

Specimen Collection - John Doe - **5555 - CCF/Specimen ID Number IS123456789**

Drug Name	Result
Amphetamines	NEGATIVE
Barbiturates	NEGATIVE
Benzodiazepines	NEGATIVE
Cocaine	NEGATIVE
Marijuana	NON-NEGATIVE
MDMA	
Methadone	
Methamphetamines	
Opiates	
Oxycodones	

Validity

Validity Name

POCT Non-Negative Results

All non-negative results are confirmed by GC/MS or LC/MS/MS. For more information, please contact the laboratory.

Send this specimen to the lab for confirmation testing

Lab/Vendor: Quest Diagnostics

CCF/Specimen ID Number: IS123456789

End test no confirmation available at this time

SUBMIT A COLLECTION EVENT

⚠ ATTENTION YOU ARE ABOUT TO SUBMIT A COLLECTION EVENT

Donor Name: Doe, John

Date of Collection: Thu Jan 11 2018 10:16:03 GMT-0700 (Mountain Standard Time)

Collector Name: KAREN GARCIA

Product Name: POCT Collection

Description: NEXSCREEN CLIA 10 PANEL

RESULT = NON-NEGATIVE

Printing CCF1 Testing Facility Copy Window

New i3ccf Order

Use this wizard to order, collect, and result an i3ccf order. **All required fields are in bold**

Collection Information

Company: ⚠ This field is required

Company Location:

Package: **Reason for Test:**

Order Expiration: Location Code: Federal Agency: Observed Collection Required:

Collection Site

Collection Site: **View/Print CCF 1: Testing Facility Copy** **View/Print CCF 2: Donor Signed Copy (2-5)** **Cancel: Print Later**

ABC Company Demo Account Collection S

Donor Information


First Name: **Middle Name:** **Last Name:**

SSN/EID: **DOB:** **Email:**

Daytime Phone: **Evening Phone:**

CCF1 Laboratory Copy Example

Note the Affix matching label section should be present, Account # should be present and check section G. Section G should start with a Quest test code rather than list a NexScreen Device. Reach out to Operations@nexscreen.com if you have questions.

FORENSIC DRUG TESTING CUSTODY AND CONTROL FORM	
<div style="border: 1px dashed gray; padding: 10px; margin: 0 auto; width: 80%;"> <p style="text-align: center;">Affix matching donor-signed consent label with barcode here Specimen ID: IS123456789</p> </div>	
Specimen ID: IS123456789	Acct: QD1234567
	
STEP 1: COMPLETED BY COLLECTOR OR EMPLOYER REPRESENTATIVE	
<p>A. Employer Name, Address, I.D. No. ABC Company Demo Account 9501 Northfield Boulevard Denver, CO 80238 PH:555-555-5555 FX: Location Code:</p>	<p>B. MRO Name, Phone and Fax No. David Nahin M.D. i3screen 9501 Northfield Blvd.Denver, CO 80238 PH:877-585-7366 FX:855-253-5666</p>
<p>C. Donor SSN or Employee I.D. No. 555555555</p>	<p>D. Donor Name Doe, John</p>
<p>E. Donor ID Verified Photo ID</p>	<p>F. Reason for Test PRE-EMPLOYMENT</p>
<p>G. Drug Tests to be Performed: 28870N NEXSCREEN 10-PANEL CONFIRMATION</p>	<p>H. Collection Site Collector</p>
<p>ABC Company Demo Account Collection Site -i3CCF Collection 9501 Northfield Boulevard</p>	<p>Collector Phone No.: 555555555</p>
<p>STEP 2: COMPLETED BY COLLECTOR (make remarks when appropriate) Collector reads specimen temperature within 4 minutes.</p>	